



The weights for your PAL lift are located in the bottom of the pallet. Remove the restraining planks with a screwdriver to gain access to the weights.

The weight block is approximately 30 pounds. Use caution when lifting and be careful not to drop the weights on your fingers. The finger holds on each end facilitate handling.

The weights leave the factory in good condition, however, they may have shifted during shipping. Inspect the weights prior to installation to be sure the acrylic paint surface is intact. If there are any scratches in the paint surface, use the touch-up paint supplied with the lift to cover the scratches. The weights are made from cast iron. If the raw surface of the cast iron is exposed to the elements, it will develop rust which could drip onto the pool deck. Use the touch-up paint to cover any areas where the raw cast iron is exposed.

Once the weights pass your inspection, place them onto the PAL frame as shown in the diagrams. Stack the weights so that the feet on the bottom of the weights nest into the indentations on the top of the weights. There are three layers of weights, eight per layer.

Once the weights are in place, attach the restraining harness as shown in the diagrams, snapping onto the snaps underneath the frame. Then assemble the rest of the PAL according to the instructions in the User's Manual.

You should inspect your weight stack at least twice a year to check for any scratches in the acrylic surface. Use the supplied touch-up paint to cover any area of raw cast iron that is showing through the surface. You can also use any acrylic spray paint for this purpose.

Please call our Customer Service Department at 800-577-4424 if we can be of further assistance.

**NAME:**  
**DIAGRAM FOR INSTALLING WEIGHTS AND RESTRAINING HARNESS**

**QTY/UNIT**

**24 WEIGHTS**

THIS INFORMATION AND DATA CONTAINED HEREIN IS CONSIDERED TO BE THE EXCLUSIVE PROPERTY OF Rehomed International. THIS DOCUMENT MAY NOT BE COPIED OR PUBLISHED OR DISCLOSED WITHOUT THE EXPRESSED WRITTEN CONSENT OF Rehomed International.

DRAWN: LC DATE: 9-18-08  
REV: LC DATE: ---  
DRAWN: JC DATE: 9-18-08

UNSPECIFIED TOLERANCE  
.XX ± 0.01  
.XXX ± 0.005

ANGLE: ±1°  
FRAC: ±1/64



Rehomed International  
522 W. Mowry Drive  
Homesstead, FL 33030  
T: 305-247-8300  
F: 305-247-8304  
E: info@poollifts.com

NAME: P/N:

REVISION #

SERVER: customers/products/weights