



RehaMed International Limited

176 Worlds End Lane, Chelsfield Park,
Orpington, Kent BR6 7SS

Tel: +44 (0)1689 851413

Fax: +44 (0)1689 857765

Mobile: +44 (0)7831 693513

E-mail: Martinsthompson@aol.com

WARRANTY INFORMATION

All Aquatic Lifting Systems have a three (3) year warranty on the frame, excluding the powder coated paint finish, which may become scratched with normal use. All electronic and motor components, with the exception of batteries, have a full two (2) year warranty. Within the warranty period, the manufacturer will repair or replace any item deemed to be found defective. Lift batteries come with a one-year pro-rated warranty. During the first 90 days of ownership, batteries will be covered 100%. If a battery failure occurs between day 91 and day 365, batteries are covered at 50% of the original cost. Normal maintenance and care of the unit, including charging of the battery when not in use is recommended. Do not store the unit, battery or components near or around chemicals.

The manufacturer shall not be responsible for the cost of removal or replacement of any defective product, nor for any other expenses or for damages which might be incurred in such removal and replacement.

This warranty specifically excludes fading of materials, microbiological staining and rust or corrosion of any metallic products or parts. Refer to the manufacturer's care and maintenance instructions for regular maintenance and cleaning of products.

This warranty relates only to defects in materials and workmanship and does not include damage or failure resulting from other causes, including, but not limited to Acts of God, misuse or abuse, accident or negligence, fire, improper assembly or installation, chipping or flaking of powder or vinyl coatings, or ice damage. Damage induced by the improper use of chemicals is not covered by this warranty. In the event that products are altered or repaired by anyone without the prior written approval by the manufacturer, all warranties are void.

The warranty is non-transferable and is subject to terms and conditions which can be viewed in full at www.srsmith.com.

Batteries

All batteries are inspected prior to shipment and should be free from defect. See warranty policy regarding battery replacement due to defect. Batteries have a normal lifespan of between 2-3 years, depending on use and care.

Battery should be charged daily. Do not allow battery to fully discharge, as it will shorten battery life. Whenever battery is removed from Control Box or lift is not in use it is important to connect battery to charger. This ensures battery is fully charged and ready for use.

Do not drop the battery, as it could cause the unit to fail. If the battery case is cracked do not use and replace. Battery has contacts located on the bottom of the housing. Do not place on a conductive surface that can cause a short. During temperature extremes beyond the range of 41 F (5 C) to 104 F (40 C) remove battery and place in a controlled environment or battery life can be shortened.

Keep the Console/Battery Cover installed at all times to prevent moisture from collecting on the Control Box and battery. Excessive moisture collection can affect battery and lift performance and could lead to battery failure and/or the lift failing to operate.

Battery Disposal

The battery is recyclable and shall be disposed of in accordance with applicable local, state/provincial or federal/national regulations.

*Manufacturers of Lifting Devices
for any Aquatic Occasion!*

PAL- Portable Aquatic Lift

SPLASH! Aquatic Lift

SPLASH! Mariner

Voyager Series

Lifts for Cruise Ships

Registered in England and Wales Number: 05195610