



BOM			
No.	P/N	DESCRIPTION	QTY
1	200-1000	PAL FRAME ASSEMBLY	1
2	150-1000	MAST ASSEMBLY	1
3	200-5000	PAL HOUSING	1
4	200-7000	HANDLE	1
5	100-1000	CONTROL BOX	1
6	100-2000	BATTERY	1
7	110-3100	WEIGHT	24
8	160-1000	SEAT ASSEMBLY	1
9	800-2039	4 ARM KNOB	1
10	800-2042	HBOLT, $\frac{5}{16}$ -16 x 2 $\frac{1}{4}$	1
11	800-2003	HNUT, NYLON, $\frac{1}{2}$ -13	2
12	800-1160	SPLIT LOCK WASHER, $\frac{1}{2}$	2
13	200-5400	PAL DECALS	3



RehaMed International  
 522 W. Mowry Dr.  
 Homestead, FL 33030  
 T: 305-247-8300  
 F: 305-247-8304  
 E: info@poollifts.com

Description: **PAL FINAL ASSEMBLY**

THIS INFORMATION AND DATA CONTAINED  
 HERIN ARE THE EXCLUSIVE PROPERTY OF  
 RehaMed International. THIS DOCUMENT MAY  
 NOT BE COPIED OR PUBLISHED OR  
 DISCLOSED WITHOUT THE EXPRESSED  
 WRITTEN CONSENT OF RehaMed International.

Drawn: SY 08-30-09  
 Checked: LC 09-11-09  
 Approved: LC 09-11-09

P/N: 200-0000  
 Scale: NTS  
 Date: 09-11-09  
 File: eng/cust/Prod-009/Explo/200-0000 PAL  
 Page: 1 of 1